Compex® WIRELESS

User Guide

1 INTRODUCTION	3
2 WARNINGS	3
Contraindications	3
Safety measures	3
3 PRESENTATION	5
Warranty	6
Maintenance	6
Storage and transportation conditions	6
Use conditions	6
Elimination	6
Standards	6
Patents	7
Standardized symbols	7
Technical characteristics	7
4 HOW DOES ELECTROSTIMULATION WORK?	8
5 USAGE GUIDELINES	9
Placement of electrodes	9
Stimulation positions	9
Adjusting stimulation energies	10
Progression in the levels	10
6 HOW IT WORKS	11
1 Description of the device	2
2 1st use: DEMO mode	3
3 Synchronize your remote	4
4 Start with your device	5
4.1 Select schedule, objective, setting or programme menu	6
4.2 Connect the modules to the electrodes	9
4.3 Prepare your session	9
4.4 Start your session	10
5 Charging	12
6 Problems and solutions	18
7 OBJECTIVE	21
8 PROGRAMMES	21

1 INTRODUCTION

Please read this manually carefully prior to using your Compex Wireless unit. You are strongly advised to read Chapter 2 of this manual. The Compex Wireless is a stimulator intended for muscular training and pain relief. Any person excepted those mentioned in the Chapter 2 "Warnings" can use the Compex Wireless device.

2 WARNINGS

No modification of this equipment is allowed!

CONTRAINDICATIONS

- Cardiac stimulator (pacemaker)
- Epilepsy
- Pregnancy (do not use on abdominal region)
- Serious arterial circulation disorders in lower limbs
- Abdominal or inguinal hernia

OSTEOSYNTHESIS EQUIPMENT

The presence of osteosynthesis equipment (metallic equipment in contact with the bone: pins, screws, plates, prostheses, etc.) is not a contraindication. The electrical currents of the Compex Wireless device are specially designed to have no harmful effect on osteosynthesis equipment.

Important!

- Do not use the Compex Wireless programmes if you have sensitivity problems.
- Never use the Compex Wireless unit for prolonged periods without medical advice.
- Consult your doctor if you want to use the device for rehabilitation purposes.
- People unable to express themselves cannot use the device.
- Consult your doctor if you are in any doubt whatsoever.
- Read this manual carefully, in particular chapter 7, which provide information concerning the effects and indications of each stimulation programme

SAFETY MEASURES

- Do not apply stimulation near the area of an implant, such as cochlear, pacemakers, skeletal or electrical.
- Do not apply stimulation in the vicinity of metal. Remove jewellery, body piercings, buckles or any other removable metallic product or device in the area of stimulation.
- If pregnant or menstruating, do not place electrodes directly over the uterus or connect pairs of electrodes across the abdomen.
- Do not use the Compex Wireless device in water or in a humid atmosphere (sauna, hydrotherapy, etc.).
- Do not use the Compex Wireless device in an oxygen-rich atmosphere.
- Never carry out an initial stimulation session on a person who is standing. The first five minutes of stimulation must always be performed on a person who is sitting or lying down. In rare instances, very highly strung people may have a vagal reaction. This is of psychological origin and is connected with a fear of the muscle stimulation as well as surprise at seeing one of their muscle contracts without having intentionally contracted it themselves. This vagal reaction produces a feeling of weakness with a fainting tendency, slowing down of the heart rate and a reduction in arterial blood pressure. If this does occur, all that is required is to stop the stimulation and lie down with the legs raised for the time it takes for the feeling of weakness to disappear (5 to 10 minutes).
- Never allow movement resulting from muscular contraction during a stimulation session. You should always stimulate isometrically; this means that the extremities of the limb in which a muscle is being stimulated must be firmly fixed, so as to prevent the movement that results from contraction.
- Do not use the Compex Wireless device if you are connected to a high-frequency surgical instrument as this could cause skin irritation or burns under the electrodes.
- Do not use the Compex Wireless unit within one metre of short wave or microwave devices as this could alter the currents generated by the stimulator. If you are in any doubt as to the use of the stimulator close to another medical device, seek advice from the manufacturer of the latter or from your doctor.
- Do not use the Compex Wireless device in areas where unprotected devices are used to emit electromagnetic radiation. Portable communications equipment can interfere with electrical medical equipment.
- Do not disconnect any module during a session. Switch it off first.
- Always use the AC adaptor provided by Compex to recharge the stimulator.
- Never use the Compex Wireless unit or the AC adaptor if it is damaged or if the enclosures are open. There is a risk of electric discharge.
- Disconnect the AC adaptor immediately if there is abnormal heating or smell, or if smoke comes from the AC adaptor or the unit.
- Do not place the docking station in a confined space for recharging (carrying case, drawer, etc.). There is a risk of fire or electric shock.
- Keep the Compex Wireless unit and its accessories out of reach of children.
- Do not allow any foreign bodies (soil, water, metal, etc.) to penetrate the device and the AC adaptor.
- Do not spread any liquid on the unit.
- Sudden temperature changes can cause condensation to build up inside the stimulator. Use the device only once it has reached ambient temperature.
- Do not use the stimulator while driving or operating machinery.
- Do not apply stimulation while sleeping.
- Do not use the Compex Wireless unit at altitudes of over 3,000 metres.

PRECAUTIONS WHEN USING THE COMPLEX WIRELESS UNIT

- After trauma or a recent operation (less than 6 months)
- Muscular atrophy
- Persistent pain
- Need for muscular rehabilitation

WHERE NEVER APPLY THE ELECTRODES

- In the vicinity of the head.
- At the front and side areas of the neck.
- Contralaterally, i.e. do not use two poles connected to the same channel on opposite sides of the body.
- On or near skin lesions of any kind (wounds, swelling, burns, irritation, eczema, cancer lesion, etc.).

PRECAUTIONS WHEN USING THE ELECTRODES

- Only use electrodes supplied by Compex. Other electrodes may have electrical properties that are unsuitable for the Compex Wireless unit. Caution should always be exercised with current densities greater than 2mA/cm².
- Always turn off the stimulator before moving or removing any electrode during a session.
- Do not place the electrodes in water.
- Do not apply solvents of any kind to the electrodes.
- Do not attempt to place electrodes on any part of the body not directly visible without assistance.
- For best results, wash and clean the skin of any oil and dry it before attaching the electrodes.
- Attach the electrodes in such a way that their entire surface is in contact with the skin.
- For obvious reasons of hygiene, each user must have his/her own electrodes set. Do not use the same electrodes on different people.
- Never use a set of adhesive electrodes for more than 15 sessions as the quality of the contact between the electrode and the skin, which is important for the comfort and effectiveness of stimulation, progressively degrades.
- Some people with very sensitive skin may experience redness under the electrodes after a session. Generally, this redness is completely harmless and disappears after 10 to 20 minutes. However never start another stimulation session in the same area if the redness is still visible.
- Please consult electrodes use and storage instructions data on electrodes packaging.

3 PRESENTATION

RECEPTION OF EQUIPMENT AND ACCESSORIES

Your kit contains:

- A. 1 remote control (6522026)
- B. 4 modules (984337)
- C. 1 docking station (683121)
- D. 1 AC adaptor (693018)
- E. 1 USB cable (601163)
- F. 2 bags of small electrodes (5x5 cm) (42215)
- G. 2 bags of large electrodes (5x10 cm 2snap connections) (42216)
- H. 2 bags of large electrodes (5x10 cm 1snap connection) (42222)
- I. 1 quick start guide (M321318)
- J. 1 travel pouch (680029)
- K. 1 Warning leaflet (xxxxxx)

WARRANTY

See the enclosed leaflet.

MAINTENANCE

Clean with a soft cloth and an alcohol based, solvent-free cleaning product. Use only a minimum amount of liquid when cleaning the unit. Do not disassemble the remote control, modules, docking station or AC adaptor as all these contain high voltage components that could cause electric shock. This may only be done by persons or repair services authorized by Compex. Your device does not require calibration. If your device contains parts that seem worn or defective, please contact your nearest Compex customer service centre.

STORAGE AND TRANSPORTATION CONDITIONS

The unit must be stored and transported in accordance with the following conditions:

Temperature: - 20° C to 45° C
Maximum relative humidity: 75 %
Atmospheric pressure: 700 hPa to 1060 hPa

USE CONDITIONS

Temperature: 0° C to 40° C
Relative humidity: 30 % to 75 %
Atmospheric pressure: 700 hPa to 1060 hPa
Do not use in areas where there is a risk of explosion.

ELIMINATION

Battery must be disposed in compliance with relevant national regulatory requirements. Any products bearing the WEEE marking (a crossed-out wheeled bin) must be separated from ordinary household waste and sent to special collection facilities for recycling and recovery.

STANDARDS

To guarantee your safety, the device has been designed and manufactured in compliance with the requirements of the European Medical Device Directive 93/42/EEC as amended.

The Compex Wireless unit also complies with the general safety standard for electromedical equipments, EN 60601- 1. It also complies with the electromagnetic compatibility standard, EN 60601-1-2 and the standard on particular safety requirements for nerve and muscle stimulators, EN 60601-2-10.

Current international standards require that a warning be given concerning the application of electrodes to the thorax (increased risk of cardiac fibrillation).

The device also complies with Directive 2002/96/EEC on waste electrical and electronic equipment (WEEE).

PATENTS

The Compex Wireless device incorporates several innovations with patents pending or delivered.

STANDARDIZED SYMBOLS


Caution / Follow the user manual or operating instructions.

The Compex Wireless unit is a class II device with internal electric power and type BF applied parts.

Name and address of the manufacturer and manufacturing date.

Name and address of the authorised representative in the European Community.

This device must be separated from ordinary household waste and sent to special collection facilities for recycling and recovery.

The On/Off button is a multi-function button.

Non-ionizing radiation

Keep away from sunlight

Keep dry

Latex Free

TECHNICAL CHARACTERISTICS

A General information

Remote control battery: Lithium Polymer (LiPo) rechargeable 3.7[V] / ≥ 1500 [mAh].

Module battery: Lithium Polymer (LiPo) rechargeable 3.7[V] / ≥ 450 [mAh].

AC adaptor: Only AC adaptor 5[V] / 3.5[A] bearing the part number 6490xx can be used to recharge the Compex Wireless unit.

B Neurostimulation

All electrical specifications are given for an impedance of 500-1000 ohms per channel.

Channels: four independent and individually adjustable channels that are electrically isolated from each other.

Impulse shape: constant rectangular current with pulse compensation to eliminate any direct current component to prevent residual polarisation at skin level.

Maximum pulse intensity: 120 mA.

Pulse intensity increments: manual adjustment of stimulation intensity from 0 to 999 (energy) in minimum increments of 0.25 mA.

Pulse width: 50 to 400 μ s.

Maximum electrical charge per pulse: 96 microcoulombs ($2 \times 48 \mu$ C, compensated).

Standard pulse ramp-up time: 3 μ s (20%-80% of maximum current).

Pulse frequency: 1 to 150 Hz.

C RF data

Frequency band of transmission: 2.4[GHz] ISM

The type and frequency characteristics of the modulation: GFSK, +/-320[kHz] deviation

The effective radiated power: TBD [dBm] Waiting for IEC result

D Information related to electromagnetic compatibility (EMC)

The Compex Wireless unit is designed to be used in typical domestic environments and is approved according to the EMC safety standard of EN 60601-1-2.

This device emits very low levels in the radio frequency (RF) interval and is therefore not likely to cause any interference with nearby electronic equipment (radios, computers telephones etc.).

The Compex Wireless unit is designed to withstand foreseeable disturbances originating from electrostatic discharges, mains supply magnetic fields or radio frequency transmitters.

Despite this, it is not possible to guarantee that the stimulator is not affected by strong RF (radio frequency) fields emitting from others sources.

For more detailed information regarding electromagnetic emission and immunity, please contact Compex.

4 HOW DOES ELECTROSTIMULATION WORK?

The principle of electrostimulation is to stimulate nerve fibres by means of electrical impulses transmitted by electrodes. The electrical pulses generated by stimulators designed by Compex are high quality pulses - offering safety, comfort and efficiency - which can stimulate different types of nerve fibres:

1. The motor nerves, to stimulate a muscular response. The quantity and the benefits obtained depend on the stimulation parameters and this is known as electromuscular stimulation (EMS).
2. Certain types of sensitive nerve fibres to obtain analgesic or pain-relieving effects.

1. Motor nerve stimulation (EMS)

In voluntary activity, the order for muscular work comes from the brain, which sends a command to the nerve fibres in the form of an electrical signal. This signal is then transmitted to the muscle fibres, which contract. The principle of electrostimulation accurately reproduces the process observed during a voluntary contraction. The stimulator sends an electrical current impulse to the nerve fibres, exciting them. This excitation is then transmitted to the muscle fibres causing a basic mechanical response (= muscular twitch). The latter constitutes the basic requirement for muscular contraction. This muscular response is completely identical to muscular work controlled by the brain. In other words, the muscle cannot distinguish whether the command comes from the brain or from the stimulator. The parameters of the programmes (number of impulses per second, contraction time, rest time, total programme time) subject the muscles to different types of work, depending on the muscle fibres. In fact, different types of muscle fibres may be distinguished according to their respective contraction speed: slow, intermediate and fast fibres. The fast fibres will obviously predominate in a sprinter, while a marathon runner will have more slow fibres. With a good knowledge of human physiology and a perfect mastery of the stimulation parameters of the various programmes, muscular work can be directed very precisely towards the desired goal (muscular reinforcement, increased blood flow, firming up, etc.).

2. Stimulation of the sensory nerves

The electrical impulses can also excite the sensory nerve fibres to obtain an analgesic or pain relieving effect. The stimulation of the tactile sensory nerve fibres blocks the transmission of pain by the nervous system. The stimulation of another type of sensory fibres creates an increase in the production of endorphins and, therefore, a reduction in pain. With pain relief programmes, electrostimulation can be used to treat localised acute or chronic pain as well as muscle pain.

Attention : Do not use the pain relief programmes for a long period without medical advice.

BENEFITS OF ELECTROSTIMULATION

Electrostimulation is a very effective way to make your muscles work:

- With significant improvement of different muscle qualities
- Without cardio-vascular or mental fatigue
- With limited stress on the joints and tendons. Electrostimulation thus allows a greater quantity of work to be done by the muscles compared with voluntary activity.

To be effective, this work must involve the greatest possible number of muscle fibres. The number of fibres that are working depends on the stimulation energy. It will therefore be necessary to use the maximum tolerable energy. The user controls this aspect of stimulation. The higher the stimulation energy, the greater the number of muscle fibres that are working and, therefore, the more significant the progress achieved. To maximize your results, Compex recommends that you complement your electrostimulation sessions with other efforts, like:

- Regular physical exercise
- Proper and healthy nutrition
- Balanced lifestyle

5 USAGE GUIDELINES

The usage guidelines presented in this section should be considered as general rules.

PLACEMENT OF ELECTRODES

For optimal results, use the recommended electrode placement. A different electrode must be connected to each module.

Note: It is possible and normal to have an electrode arrangement that leaves one electrode connection free.

For programmes which force muscles to contract, it is important to place the electrode attached to the pod with the on/off button on the motor point of the muscle. It is important to choose the right size of electrodes (large or small) and correctly position these on the muscle group you want to stimulate so as to ensure the efficacy of the treatment. Therefore, always use the size of electrodes shown in the pictures. Unless you have other specific medical instructions, always follow the placement directions in the pictures. Where necessary, look for the best possible position by moving the positive electrode over the muscle until you find the point that will produce the best contraction or the most comfort for you.

Note: Approximate placing of the electrodes makes for a less efficacious session, but is still risk-free.

Compex disclaims all responsibility for consequences arising from electrodes placed in other positions.

Electrode placement can be found on the www.compexwireless.com/coach on the Objective or a Programme selector or directly in the remote control under the  icon when you select an Objective or a Programme.

STIMULATION POSITIONS

This position will vary depending on the position of the electrodes, the muscle group you wish to stimulate, and the programme you are using. For programmes which involve powerful muscular contractions, the muscle should always be stimulated in an isometric way. You must therefore fix the extremities of your limbs securely. In this way, you provide maximum resistance to the movement and prevent any shortening of the muscle during the contraction, which could create cramp pains and serious stiffness after the session. For example, when stimulating the quadriceps, the user should be in a seated position with the ankles fixed with straps to prevent extension of the knees.

ADJUSTING STIMULATION ENERGIES

In a stimulated muscle, the number of recruited fibres depends on the stimulation energy. For programmes which involve powerful muscular contractions, you must therefore use maximum stimulation energies (up to 999), always at the limit that you can endure, in order to recruit the maximum number of fibres.

For programmes which involve only muscle twitches, you should progressively increase stimulation energies until you obtain clearly visible muscular twitches.

PROGRESSION IN THE LEVELS

In general, it is not advisable to go through the different levels quickly with the intention of reaching level 3 as fast as possible. In fact, the different levels correspond to progress with electrostimulation. The simplest and most usual procedure is to start with level 1 and raise the level when changing to a new stimulation cycle. At the end of a cycle, you may either start a new cycle at the next level up or do some maintenance at the rate of 1 session

6 HOW IT WORKS

You are strongly advised to read carefully the contraindications and safety measures described at the beginning of this manual in Chapter 2 prior to using your unit.

2 DESCRIPTION OF THE DEVICE

REMOTE CONTROL

- A** On/Off button
- B** 4 multi-functions buttons:
 - Function related to the pictograms that are on the screen (i.e. info, back, home, etc.)
 - Selection of the desired stimulation channel to increase or decrease the energy
- C** Multi-directions (up-down-left-right) pad to select items
- D** Validation button
- E** Socket for the USB cable or the docking connector

MODULE

- A** On/Off button (press briefly to switch on but press and hold to switch off)
 - LED blinking green: ready to use
 - LED blinking yellow: on stimulation
- B** Groove to wind up the cable

DOCKING STATION

- A** Connector for the remote control
- B** Slit to open the docking top
- C** Slot to position the module to recharge it
- D** Socket for the AC adaptor

CONNECTIONS

Connect the AC adaptor supplied with your device to the docking station, then plug the AC adaptor into a socket. We strongly recommend you to fully charge the battery of the remote control and modules before using them for the first time as this will improve its performance and life span (see chapter Recharging).

2 1ST USE: DEMO MODE

SELECT THE LANGUAGE

DEMO mode is available as long as you don't pair your device with your computer. DEMO mode gives access to several programmes only and not to all the possibilities of the Compex Wireless (empty Calendar and no Objective).

Compex advise you to go to the www.compexwireless.com/coach and create your own account in order to select preset objective and to follow your progression on your schedule or to download new programmes.

DEMO MODE

Select your language
Validate your choice

demo Press DEMO button

3 SYNCHRONIZE YOUR REMOTE

CREATE YOUR PERSONAL ACCOUNT

If you want to benefit from all the possibilities of the product (download training planning, create your own objective, follow your progression, add new programmes, etc.) you should first create your personal account on www.compexwireless.com/coach website. Follow the instruction on the website.

If you don't want to create your personal account you can use the product in DEMO mode (see point 1st use: DEMO mode)

PAIR YOUR REMOTE CONTROL WITH YOUR ACCOUNT

As soon as your account is created you can log in your account and pair your remote control.

To do so please follow the instruction on the website www.compexwireless.com/coach

- A** Download and install the software on your computer
- B** Click on the "Pair your remote" button and type the serial nb of the remote control
- C** Connect the remote control and press the round button of the channel 4 (the one on the right of the remote control) 
- D** Type on the website the code shows on your remote control

4 START WITH YOUR DEVICE

As soon as you have synchronized your remote control with your account and switch on your device, the first screen shows the today task. Probably you will have a day off but if you have already entered an Objective you will see some tasks. From there you can also start a session (see point Select an objective or a programme)


Info Menu, press and hold

Select your language

Validate your choice

Previous screen

4.1 SELECT SCHEDULE, OBJECTIVE, SETTING OR PROGRAM MENU

Select your item

Validate your choice

1. SCHEDULE

Here you can find your schedule with the session you did or you will have to do.

- Green rectangle below the day means that you have done what is planned
- Orange rectangle below the day means that you haven't done what was planned in the past
- Grey rectangle below the day means that you have planned something in the future

Validate your choice

Select your item

Back to main menu

Validate your choice

Select your language

Previous screen

Info Menu, press and hold

By selecting a day you have the details of this day (Objective of Programme planned) and from there you can also start a session (see point Select an objective or a programme)

2. OBJECTIVE

Here you find the objective you have selected on the website and downloaded on your device.
See point Select an objective or a programme on how to navigate in the menu

3. SETTINGS

On the settings menu you can set several parameters like Backlight, Volume, mi-Scan option, Language, etc.

Validate your choice

Select your item

Back to main menu

4. PROGRAMMES

The programmes menu includes the list of your programmes.

In the DEMO mode only several programmes have been installed. As soon as you pair your remote control with your web account you will be able to change this list (add or remove programmes).
See point Select an objective or a programme on how to navigate in the menu.

SELECT AN OBJECTIVE OR A PROGRAMME

The selection of an objective or a programme is very similar.

SELECT AN OBJECTIVE

Validate your choice
Select your item
Back to main menu

SELECT A PROGRAM

Validate your choice
Select your item
Back to main menu
Info menu

WHEN AN OBJECTIVE HAS BEEN SELECTED YOU CAN SEE THE DETAILS OF ALL THE SESSIONS THAT ARE PART OF THIS OBJECTIVE

Validate your choice
Select your item
Back to the previous screen

See the electrodes placement and body position related to your objective. For programmes you have all the possible electrode placements. See them with the up/down multi-direction pad.

At this stage you can place the electrodes on your body and connect the modules to the electrodes.

Info menu, keep it pressed
N.B.: For voluntary session (cycling in this example) you have to validate yourself the training by pressing the button under this icon

* 1. Body position
2. Electrode placement

4.2 CONNECT THE MODULES TO THE ELECTRODES

Stick the electrodes on your skin. The fixation of the modules on the electrodes is made from the side. Slide the module on the electrode snap until you hear a click

To remove the modules from the electrodes simply do it on the other way.

4.3 PREPARE YOUR SESSION

WARMUP AND COOL DOWN

This screen is not available for all the programmes.

You can choose if you want the warmup phase before starting the workout session and if you want the cool down at the end of the session.
N.B.: Compex advise you to keep the warmup unless you have already performed a voluntary warmup and the cool down phase.

Confirm or not the warm up or cool down phase
Previous screen
Validate your choice

4.4 CONNECT ELECTRODE TO MODULE

SWITCH ON THE MODULE BEFORE STARTING THE STIMULATION


Just before starting the stimulation and the mi-scan measure you have to switch on one module after the other. As soon as the 1st module is detected by the remote control, you have the possibility to switch on the 2nd module and so on until the 4th module. If you want to use only 2 modules for your session just press the START button after the detection of the second module.

BEGINNING OF THE STIMULATION SESSION: MI-SCAN MEASURE

This function adapts the electrostimulation session to the physiology of each user. Just before starting the work session, mi-SCAN tests the muscle where you put the electrodes and automatically adjusts the stimulator parameters to the excitability of this body area. It is a truly personalised measurement. And in order to get the best efficacy and comfort of your session, you are advised to perform the mi-SCAN measure every time. This function is implemented at the start of the programme by a short sequence during which measurements are taken. During this phase you must remain still and relaxed. When the test is over you can start the stimulation session.

DISPLAY EXPLANATION

DURING STIMULATION SESSION

The stimulation starts at 0:00.

N.B.: Active rest phase stimulation energies are automatically set at 50% of contraction energies. These can be modified during the rest phase. Once modified, they will be totally independent of the contraction energies.

PAUSE MODE

N.B.: The session will resume at 80% of the energy levels that were in use prior to the interruption.

-  Re-start the stimulation session
-  Back to Menu before Stimulation

END OF SESSION


By switching any button at the end of the session you go back on the main menu. If you want to switch your remote control off just press and hold the On/Off button.

5 CHARGING

BATTERY LEVEL

You can see the battery level of the module when you synchronize it (before the stimulation). The battery level of the remote is always indicated on the top right corner. The small green patch shows how many modules are connected to the remote.

- Module battery level
- Number of modules connected to the remote
- Remote battery level

CONNECT THE DOCKING STATION

Connect the AC adaptor supplied with your device to the docking station, then plug the AC adaptor into a socket. We strongly recommend you to fully charge the battery of the remote control and modules before using them for the first time as this will improve its performance and life span.

CHARGING THE MODULES AND THE REMOTE CONTROL


At the end of your session you are advised to store the module and the remote control in the docking in order to recharge them. Place the remote control on the connector of the docking station.

THIS SCREEN APPEARS ON THE REMOTE

Then place the module in the docking station

To do so, place the pod without the On/Off button (the **green** one) on the slot indicated in **green** and the other pod on the **blue** slot. Do the same for the others modules.

The pod without the On/Off button (the **green** one) has to be inserted in the small connector. A magnet should help you in positioning the pod as well as the groove of the pod. When it is placed correctly you should hear a click.

WHEN YOU CONNECT A DISCHARGED MODULE ON THE DOCKING A RED CIRCLE APPEARS. AS LONG AS THE RED CIRCLE IS ON, THE MODULE IS CHARGING. AS SOON AS IT DISAPPEARS THAT MEANS THE MODULE IS CHARGED

N.B.: If the module is charged when you place it on the docking no red circle will appear.

As soon as the modules and the remote control are fully charged the actual time is displayed.

6 PROBLEMS AND SOLUTIONS

ELECTRODE FAULT

The remote control displays the symbol of an electrode disconnected on the channel where the problem has been detected (in the example above, error in channel 1).

- Check that electrodes are correctly connected to the module.
- Check whether the electrodes are old, worn, and/or the contact is poor: try using new electrodes.

MODULE OUT OF RANGE

The remote control displays the symbol out of range on the channel where the problem has been detected (in the example above, error in channel 1).

- Check that module and remote control are close to less than 1 meter.
- Be sure that you are not in a remote area without any obstacle that can reflect the signal.
- Be sure that you are in a place where the signal of the remote control can be reflected

SYNCHRONISATION ISSUE

If the synchronisation process has been interrupted or if for another reason you cannot synchronise (remote disconnected, power failure, etc.), the remote control can display but not always the screen above.

- Simply reconnect the remote control to the PC and relaunch the synchronisation process by pushing the round button of the channel 4 (the one on the right of the remote control) 

MODULE RUNS OUT OF BATTERY

During the stimulation a module can run out of battery. In this case the symbol of an empty battery is displayed on the channel where the problem has been detected (in the example above, in channel 4).

- You can stop the stimulation and recharge the empty module.
- You can drop the empty module off and continue your session without it.

MODULE LED BEHAVIOUR

LED blinking green and red alternatively: the module is not in the range of the remote control.

- Be sure that the remote control is switched on.
- Check that module and remote control are close to less than 1 meter.

LED IS ALWAYS RED

- Check that module is charged.
- Try to reboot the remote control and the modules.
- If despite this the LED is still red, please contact the consumer service approved by Compex.

STIMULATION IS NOT PRODUCING THE USUAL SENSATION

- Check that all the settings are correct and ensure the electrodes are positioned properly.
- Change the positioning of the electrodes slightly.

THE STIMULATION EFFECT CAUSES DISCOMFORT

- The electrodes are beginning to lose adhesion and no longer provide good contact on the skin.
- The electrodes are worn and need to be replaced.
- Slightly change the positioning of the electrodes.

THE UNIT IS NOT WORKING

- Check that the modules and the remote control are charged.
- Try to reboot the remote control and the modules.
- If despite this the unit is still not working, please contact the consumer service approved by Compex.

7 OBJECTIVE

With the DEMO mode there is no Objective available in the Compex Wireless. To download Objective you should go to www.compexwireless.com/coach.

Example of Objective available:

- Prepare a marathon
- Prepare a half marathon
- Prepare a cycling race
- Increase your strength
- Treat your low back pain
- Prevent your ankle twist
- ...

8 PROGRAMMES

The Compex Wireless is set up with the following programmes. As soon as you create an account on www.compexwireless.com/coach and synchronise your remote with this account you will be able to download new programmes.

Programmes that induce muscular contraction like Endurance, Resistance, Strength, etc. are proposed with 3 different levels. To better understand how to use the levels see chapter V. Usage guidelines.

PROGRAMMES	EFFECTS	USES
Capillarization	Very marked increase in blood flow Development of capillaries	In the period before competition, for those participating in endurance or resistance sports To improve endurance in those whose fitness level is not high
Toning massage	Activation of blood circulation Recovery of muscular contractile properties Invigorating effect	To prepare the muscles in an ideal manner before an unusual/one-time physical activity
Endurance	Improved absorption of oxygen by the stimulated muscles Improved performance for endurance sports	For athletes who wish to improve their performance in long-duration sporting events
Resistance	Improved lactic capacity of muscles Improved performance for conscious resistance sports	For athletes who wish to increase their capacity to obtain intense and prolonged effort
Strength	Increased maximum strength Increased rate of muscular contraction	For athletes practising a discipline that requires strength and speed
Explosive strength	Increases the speed at which the level of strength is attained Improved efficiency of explosive actions (jumping, sprinting, etc.)	For athletes practising a discipline in which explosive strength is an important factor in performance
Hypertrophy	Increased muscle volume Increased muscle resistance	For athletes who wish to increase their muscle mass
Fartlek	Training and preparation of muscles for all kinds of muscular work (endurance, resistance, strength, explosive strength) through different working sequences	At the beginning of the season to "re-set" the muscles after a rest period and before more intensive and specific training During the season for those who do not have any preferences for a single kind of performance and prefer to exercise their muscles with different work-out methods
Active recovery	Marked increase in blood flow Accelerates the elimination of waste from the muscular contraction Endorphinic effect Relaxing effect	To improve and accelerate the muscle recovery after an intensive exercise To be used during the first 3 hours after each session of intensive training or after a competition
Relaxing massage	Decreased muscular tension Drainage of the toxins responsible for the exaggerated increase of muscular tone Effect of well-being and relaxation	To eliminate uncomfortable or painful sensations, following an exaggerated increase of muscular tone
Anti-stress massage	Improvement of tissue vascularisation Decreased muscular tension associated with stress Return to calm	To restore a sensation of well-being after stressful situations
Heavy legs	Increased venous return Increased supply of oxygen to muscles Drop in muscular tension Elimination of tendency to have cramp	To restore a sensation of well-being after stressful situations

PROGRAMMES	EFFECTS	USES
Modulated TENS	Blocks transmission of pain by the nervous system	To remove the feeling of heaviness in the legs arising in unusual situations (standing for long periods, heat, hormonal imbalance linked to the menstrual cycle, etc.)
Endorphinic	Analgesic action through the release of endorphins Increased blood flow	All acute or chronic localized pain
Decontracture	Reduced muscular tension Relaxing effect	To combat chronic muscular pain
Disuse atrophy	Reactivation of the trophic action of muscle fibres damaged during atrophy	To combat recent and localized muscular pain
Muscle growth	Increased diameter and capacity of muscle fibres damaged during under-activity or inactivity imposed by some kind of pathological problem Restoration of muscle volume	Treatment to combat any reduction of muscle volume: • As a result of trauma necessitating immobilisation • Accompanying degenerative lesions to joints
Reinforcement	Increase the strength of a muscle previously atrophied Increase the strength of a muscle affected by a pathological process	After a period of using the Disuse atrophy treatment, as soon as the muscle shows signs of a slight recovery of volume or tone Until the virtually complete restoration of muscle volume
Muscle building	Improved muscle trophism Balanced increase of muscular tone and volume	At the end of rehabilitation, when the muscle has regained normal muscle volume From the onset of rehabilitation for non-atrophied muscles
Toning	To tone the muscles	To be used to begin with in order to tone up and prepare the muscles before more intensive firming work
Firming	To regain muscle firmness and restore the support function of the muscles	To be used as the main treatment for muscle firming
Shaping	To define and sculpt the body when the muscles are already firm	To be used when the firming phase is finished

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